## Device Modification Analysis Summary

Ref no.: DMA 2004/09

### Part I. Product Information/Description of Change

**Product:** DRYPIX 4000 Medical Dry Laser Imager

**Predicate Device Name & 510(k) #:** DRYPIX 7000/5000 Medical Dry Laser Imagers (K033377)

#### Description of Change:

The DRYPIX 4000 represents a modification of the above cleared device. The DRYPIX 4000 uses the same laser exposure technology, thermal development technology, film conveyance layout and imaging film as the DRYPIX 7000/5000. It also has the same indications for use as the DRYPIX 7000/5000 including the printing of digital mammography images (FFDM). It differs in component layout, more accessibility for released and jammed films and it is also designed to withstand vibration and shocks in a mobile environment. Other similarities and differences are documented in the attached SE Feature Comparison table.

### Part II. Device Modification Analysis Checklist

1. Could the change or modification significantly affect the safety or effectiveness of the subject device? **No**
2. Is there a significant change to the design, material, chemical composition, energy, or manufacturing process? **No**
3. Is there a significant change to the intended use of the subject device? **No**
4. Is there a significant change to the labeling (i.e. Warnings, additional claims) which may affect the safety or effectiveness of the device? **No**

### Part III. Attached Documentation

Please list below and attach any documentation used to support this Device Modification Analysis.

1. Regulatory Action Determination Flowchart for DRYPIX 4000
2. Ruby (DryPix 4000) Specifications, dated 8/24/2004
3. Ruby (DryPix 4000) Specifications Supplement, dated 8/24/2004
4. Substantial Equivalence (SE) Feature Comparison Table

**Ref no.:** DMA 2004/09

**PLEASE CHECK ONE:**

The above information has been reviewed and it has been judged that Fuji has no pre-market obligations to the FDA for the subject device and marketing of the device may begin immediately, subject to the general controls imposed by FDA for such a device.

A more detailed narrative for this decision shall be documented either in the comment section below or in attached documentation (such as a memo, letter, etc.).

It has been determined that a pre-market obligation exists, therefore a 510(k) pre-market notification shall be filed by FMSU and cleared by FDA prior to marketing of such a device.


A more detailed narrative for this decision shall be documented either in the comment section below or in attached documentation (such as a memo, letter, etc.).

**Comments:**

The DRYPIX 4000 is considered a modification of the previously cleared DRYPIX 7000/5000 Dry Laser Imagers. The modification consists of technology and/or performance changes. Based on this evaluation, it appears that these changes (1) do not affect indications for use, (2) do not require clinical data to evaluate safety and effectiveness, and (3) do not raise new issues of safety and effectiveness.


Therefore, it is my conclusion that the DRYPIX 4000 does not change indications for use nor does it significantly affect any safety or effectiveness issues than what is already identified in the premarket notification for predicate devices, DRYPIX 7000/5000. Therefore, I have reached a conclusion that the DRYPIX 4000 represents a modification that does not require a 510(k) submittal.

**QARA SIGNATURE:**

\_\_\_\_\_  
Frank Gianelli  
Regulatory Coordinator

\_\_\_\_\_  
November 15, 2004  
Date

**Specifications of Ruby**

Classification	Customer Requirement Specification & Other Requirements		Design Input		Product Specifications (See main text for details)
	CRS No.	Requirements	Design Input No.	Items	
1. Functional Requirements	Req - Ruby - 01	Number of trays for film supply is increased after installing the main unit and two different size trays can be used simultaneously.	Req-Ruby-01-01	Number of trays for film supply	Two trays are available. The upper tier tray is normally equipped. The lower tier tray is optional (Future addition is available on the market). Any size of films can be selected according to each tray. Size of tray is set by the qualified service personnel (as Service Setting).
		A pack of films can be accumulated on the specified loading part.	Req-Ruby-01-02	Number of accumulable films on the specified loading part	Up to 150 sheets of films can be accumulated to store them.
		Printing function of annotation strings comparable with the conventional DRYPIX7000 is equipped.	Req-Ruby-01-03	Printing function of annotation strings	Annotation strings can be printed out on the upper portion and lower portion of a film. Normal character size (40 pixels for standard), twice, three times and four times each that of normal character size can be set. The annotation strings can be indicated in accordance with DICOM Annotation BOX SOP Class or by setting configuration via PC utilities. Alphanumeric, Japanese, and Chinese characters can be indicated as readable strings. Conversion function of macro characters [reception date of printing, sender's IP address, identifier (optional string corresponding to IP address), and AE-title] is employed. This system can translate to print out such annotation strings of image format in accordance with DICOM Annotation BOX SOP Class.
		The test patterns comparable with the conventional DRYPIX7000 can be output.	Req-Ruby-01-04	Test patterns printing function	The 24-steps densitometry pattern, 17-steps densitometry pattern, QC test pattern, and SMPTE test pattern can be output to control the density. The flat test pattern can be output to check the nonuniformity. The sharpness and the double density/sharpness test patterns can be output to measure the sharpness. The lattice test pattern can be output to adjust the format. The uniformity test pattern can be output to adjust the uniformity.
		Sufficient security system functions.	Req-Ruby-01-05	Security measures	The Operating System can be updated as Service function. A anti-virus software can be installed as Service function. The pattern file of the anti-virus software can be updated by users. Images stored in HDD (e.g. to reprint images and to copy any image data to the other media) can be printed out only when the registered user gives the authentication.
	Req - Ruby - 02	Sorter attaching not less than 3 bins can be mounted for future additions.	Req-Ruby-02-01	Number of sorter bins	Three-bin sorter unit is available to set as an optional for future addition. Totally, four-bin sorting together with a tray equipped with the main unit can be allowed. The numbering of the bins with the configuration setting via PC utilities is allowed in an up to down order, and vice versa.
			Req-Ruby-02-02	Number of films accumulable on each sorter bin	Up to 30 sheets of films can be accumulated. Bins should be made of suitable materials to avoid scratching the film surface when films are taken out of such bin.
			Req-Ruby-02-03	Types of Sorting	The following sorting types are available: To output films to the desired bin by designating the bin number complying with DICOM; to output films to such bin specified by each AE-title via configuration setting; and to output films to such bin by specifying each film size via configuration setting.
			Req-Ruby-02-04	When a sorter fails	When a sorter fails, the user will remove it to release films to a tray of the main unit via the control panel. When a tray of the main unit fails, films will be output to the uppermost sorter bin.
			Req-Ruby-02-05	Sorter fitting position	Fit a sorter to the film release unit.
2. Performance Requirements	Req - Ruby - 03	This System reduces restraints on installation by applying 100V power source, the space saving function, and the like more easily in comparison with DRYPIX7000.	Req - Ruby - 03 - 01	Size of Unit	Main unit only W 600×D 550×H1000 to W 630×D 600×H1100 Main unit with sorter W 600×D 550×H1300 to W 630×D 600×H1400 ※Width and depth of the floor dimension don't include a locally protruding area.
			Req - Ruby - 03 - 02	Weight	Main unit only 130 to 150Kg Main unit with sorter 135 to 180Kg
			Req - Ruby - 03 - 03	Power supply voltage	AC100 to 120V±10% 12A or less AC200 to 240V±10% 7A or less
			Req - Ruby - 03 - 04	Voltage source capacity	1.2kVA or less
			Req - Ruby - 03 - 05	Environmental working condition (operating)	Temperature 15 °C to 30 °C Humidity 15 % to 70%RH (No dew condensation) Atmospheric pressure 700 to 1060hPa
			Req - Ruby - 03 - 06	Environmental working condition (non-operating)	Temperature 0 °C to 45 °C (Avoid freezing) Humidity 10 % to 90%RH (No dew condensation) Atmospheric pressure 500 to 1060hPa
			Req - Ruby - 03 - 07	Floor (installation surface)	Floor vibration condition : Complying with Fuji's Medical Device Quality Standards ( K12N0097A ) Floor levelness : 10 mm/m in all directions ( 1/100 or less of gradient ) Floor flatness : 10mm or less
			Req - Ruby - 03 - 08	Ventilation condition	Ventilation in an area where the system is placed: 54m <sup>3</sup> /h or less per system ( Equal to or lower than that of DRYPIX7000 )
			Req - Ruby - 03 - 09	Antimagnetic field	The system normally operates at DC 100 Gauss without any problem.
	Req - Ruby - 04	This System can be installed in a Mammography Screening Mobile (When an optional component is mounted).	Req - Ruby - 04 - 01	Vibration condition when loading it on the mobile (When moving)	Random vibration: 5 to 500Hz, X, Y, and Z axis: 24 hours each Totally 72 hours (300000 km assumed) Acceleration: 1G
		Req - Ruby - 04 - 02	Vibration condition when loading it on the mobile (When operating)	References for Operating Vibration & Shock of the Medical Device Quality Standards ( K12N0097A )	

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
Oblique rear view

Panel

## Architectural Comparison between Ruby and DRYPIX 7000

- The recording synchronous heat development system is employed with the specified conveyance layout in a same manner as DRYPIX 7000.
- Up to two (2) trays can be mounted for future additions.
- Released films can be accessed from all directions (from three directions except the right hand side, when a sorter is attached).
- Jammed films are removed from the front, top or right side face.
- Control boards and power supply system shall be placed collectively to the lower part of this System.
- Optional sorter can be attached to the upper part of this System.
- Combination of adding parts and mounting a vibration proof platform can give to harmonize with the required conditions when loading this System on a mobile.

### ■ Indication for Use

Ruby is indicated for use in providing diagnostic quality medical images on film for aid in physician diagnosis, including the printing of images and associated identification information from various digital imaging source modalities, including but not limited to, Computed Tomography (CT), Magnetic Resonance Imaging (MRI), Ultrasound, Computed Radiography, Digital Radiography, Digital Mammography and Nuclear Medicine.

**Substantial Equivalence (SE) Feature Comparison  
DRYPIX 4000 to DRYPIX 7000 & DRYPIX 5000**

<b>Parameter</b>	<b>Fuji DRYPIX 4000</b>	<b>Fuji DRYPIX 7000</b>	<b>Fuji DRYPIX 5000</b>
Light source/modulation	Laser diode Direct modulation Thermal (dry) developer	Laser diode Direct modulation Thermal (dry) developer	Laser diode Direct modulation Thermal (dry) developer
Throughput	115 films/hr (14"x17") 170 films/hr (10"x14")	180 films/hr (14"x17") 240 films/hr (10"x14") 200 films/hr (8"x10")	130 films/hr (14"x17") 180 films/hr (10"x14") 180 films/hr (8"x10")
Film type	Fuji dry imaging film DI-HL (blue base), DI-HLc (clear base)	Fuji dry imaging film DI-HL (blue base), DI-HLc (clear base)	Fuji dry imaging film DI-HL (blue base), DI-HLc (clear base)
Film size	14"x17", 14"x14", 10"x14", 10"x12", 8"x10"	14"x17", 10"x14", 8"x10"	14"x17", 10"x14", 8"x10"
Film holding capacity	150 sheets + 150 sheets (2 <sup>nd</sup> film tray optional)	100 sheets + 100 sheets (2 <sup>nd</sup> , 3 <sup>rd</sup> film trays optional)	100 sheets + 100 sheets (2 <sup>nd</sup> , 3 <sup>rd</sup> film trays optional)
Film loading	Daylight loading (no magazine)	Daylight loading (no magazine)	Daylight loading (no magazine)
Sorter	Yes (optional)	Yes (optional)	Yes (optional)
Sampling raster	100 µm x 100 µm 50 µm x 50 µm	100 µm x 100 µm 50 µm x 50 µm	100 µm x 100 µm 50 µm x 50 µm
Rasters of main and auxiliary sampling	10 pixels/mm 20 pixels/mm (50 µm mode)	10 pixels/mm 20 pixels/mm (50 µm mode)	10 pixels/mm 20 pixels/mm (50 µm mode)
Maximum number of recording pixels	3520 x 4280 5080 x 7200 (50 µm mode)	3520 x 4280 5080 x 7200 (50 µm mode)	3520 x 4280 5080 x 7200 (50 µm mode)
Time to output 1 <sup>st</sup> film	80 sec	65 sec	65 sec
Effective Gray Scale	14 bits	14 bits	14 bits
Maximum Density (Dmax)	3.0 to 3.6 (configurable)	3.0 to 3.6 (configurable)	3.0 to 3.6 (configurable)
Film density control	Yes (Auto)	Yes (Auto)	Yes (Auto)
Gradation curve variability	8 independent types selectable	8 independent types selectable	8 independent types selectable
Negative/positive selection	Yes	Yes	Yes
Black/clear border changeover	Yes	Yes	Yes
Format	1, 2, 3, 4, 6, 8, 9, 12, 15, 16, 20, 24 or mixed 4, 25, 30, 35, 36,42, 48, 49, 54, 56, 60, 63, 64	1, 2, 3, 4, 6, 8, 9, 12, 15, 16, 20, 24 or mixed 4, 25, 30, 35, 36,42, 48, 49, 54, 56, 60, 63, 64	1, 2, 3, 4, 6, 8, 9, 12, 15, 16, 20, 24 or mixed 4, 25, 30, 35, 36,42, 48, 49, 54, 56, 60, 63, 64
DICOM capability	Yes	Yes	Yes
Size ("w" x "d" x "h")	630mm x 600mm x 1100mm	735mm x 680mm x 1240mm	735mm x 680mm x 1240mm
Weight	150 kg (330 lbs)	200 kg (440 lbs)	203 kg (448 lbs)
Power	100-120 Vac, 12A, 50/60Hz 200 to 240 Vac, 7A, 50/60Hz	200/220/230/240 Vac, 11/10/9.5/9A, 50/60Hz	100/110/120 Vac, 12A 50/60Hz