


FUJI MEDICAL DRY IMAGING FILM DI-HT


Fuji Medical Dry Imaging Film DI-HT is designed for use with FUJI MEDICAL DRY IMAGER DRYPIX 2000.

■ SENSITIVITY & CONTRAST

Sensitivity and contrast of the DI-HT are suitably designed for thermal-recording with FUJI MEDICAL DRY IMAGER DRYPIX 2000.




■ FILM STRUCTURE



■ AVAILABLE SIZE & QUANTITY

Film Size	Quantity / Package
① 26 × 36	100 sheets + 1 cleaning sheet
② 35 × 43	100 sheets + 1 cleaning sheet
③ 20 × 25	100 sheets + 1 cleaning sheet

■ DIMENSIONS



■ STORAGE CAUTIONS

- Make sure to keep unrecorded film in its moisture-proof bag, stored in a cool, dark place(25°C/77°F or below).
- Observe the following points in handling and storing recorded film.
Store recorded film in a cool, dry place of low temperature and low humidity. The higher the temperature and humidity, the more the density of recorded images will increase. Long-term storage at high temperatures, high humidity and/ or daylight room conditions may cause discoloration.

Results of accelerated testing indicate the following life spans for recorded film, in terms of how long it takes for a 10% rise in density to occur in the places in which image density is 1.2 at time of output.

- ① When stored at 25°C, 60%RH, approx. 30 years.
- ② When stored at 30°C, 60%RH, approx. 5 years.

If storage temperature exceeds these levels, lifespan will be as follows.

- ③ When stored at 35°C, 60%RH, approx. 6 months.
- ④ When stored at 40°C, 60%RH, approx. 1 week.

Thus, make sure not to leave recorded film in places exposed to direct sunlight, next to heating equipment, or inside automobiles.

Under the following environmental conditions, the lifespan will be approx. 5 years.

Whithin a one-year period:

- When stored 4 months at 35°C, 60%RH for 8 hours a day and 28°C, 60%RH for 16 hours a day.
- When stored 8 months at 28°C, 60%RH, 24 hours a day.

*For other details on handling of DI-HT film, refer to the operation manual.

SYMBOLS AND ABBREVIATIONS



Batch code



Expiry (Use by)

NIF

Non Interleaved Film



Store film at 10 to 25°C, at 30 to 60%RH



Store film properly shielded from X-rays, gamma rays or other penetrating radiations and the direct sun.

FUJIFILM
FUJIFILM Corporation

26-30, Nishiazabu 2-chome, Minato-ku, Tokyo 106-8620, Japan



European Authorized Representative:
FUJIFILM Europe GmbH
Heesenstrasse 31, D-40549 Düsseldorf Germany