

# DELWORKS **FIT**

## Tablet Workstation

Maximize portability and enhance efficiency

### Features

Light-weight, 10.1" Multi-touch Tablet

Designed to withstand the elements of any radiology  
(or non-radiology) environment

Up to 10 hours of battery life



# Specifications

Operating System	Windows 7 Pro 64-bit
Processor	2.0GHz, Intel Core i5
Storage/Memory	8GB RAM, 256GB SSD HD
Display	10.1" , 1920 x 1200 with LED backlighting
Keyboard/Input	10-point multi-touch supports stylus pen
Battery	Li-Ion, 10.8V 4400mAh
Dimensions	10.6" x 7.4" x 0.8" (L x W x H)
Weight	2.4 lbs.
Durability/Hazardous Location Certifications	ANSI / ISA 2.12.01-2013 MIL-STD-810G (4' drop, shock, vibration, rain, dust, sand, altitude, freeze/thaw, high/low temperature, temperature shock, humidity, explosive atmosphere)
WiFi	802.11a/b/g/n/ac



08.2015 #9000-DW-FITWorkstation



[www.delmedical.com](http://www.delmedical.com)



**New York**  
28 Calvert Street  
Harrison, NY 10528  
(Tel) 800.261.9808  
(Fax) 914.835.6111

**Chicago**  
241 Covington Drive  
Bloomington, IL 60108  
(Tel) 800.800.6006  
(Fax) 847.288.7011

**DEL MEDICAL**  
A **UMG** COMPANY