

Introducing:

evoView

DICOM Viewing Software

The ideal choice for any Chiropractic, Veterinary, Podiatry, or any type of Medical clinic looking for a diagnostic viewing and storage solution with PACS like functionality for a fraction of the cost!

Features:

- ❖ DICOM Storage SCU / SCP
- ❖ DICOM Query / Retrieve SCU
- ❖ DICOM Print SCU
- ❖ DICOM CD Burning w/ embedded DICOM viewer
- ❖ Stack / Scout view
- ❖ Thumb nail view
- ❖ Convert DICOM to JPG or BMP images
- ❖ Hanging Protocols
- ❖ Audit Trails
- ❖ Image / Study Reporting
- ❖ E-Mail Enabled

WORKLIST - [evoView]

Database Query & Retrieve DICOM CD-R Open CD Backup

* Patient ID * Sex * Patient Name * Modality * Patient Age * Study Desc
 * Report * Verify * Doctor 5/15/2008 5/15/2008

Local Backup CD Search

Sex	PID	Patient Name	Modality	Disc	No.	Series	#Image	Modality	BodyPart
F	ID2940	PATIENT2940	CR		1	1	2	CR	CHEST
M	ID2946	PATIENT2946	CR						
M	ID2930	PATIENT2930	CR						
M	ID2934	PATIENT2934	CR						
M	ID2937	PATIENT2937	CR						
M	ID2945	PATIENT2945	CR						
F	ID2935	PATIENT2935	CR						
F	ID2943	PATIENT2943	CR						
F	MS-MR-111113	Merical Standard	MR						

No.	Series	#Image	Modality	BodyPart	Study Date	Series Desc
1	1	2	CR	CHEST	2003/11/16 10:12:47	

EXAM# [60] SER ES# [151] IMAGE# [465]

start

MSPPW - [evoView]

File Edit View Image Process Annotation Tool W/L Preset LJT Preset Sharpen Preset Option Hanging Preset: Macro

1 ID2940 PATIENT2940 F 67Y 2005/04/12 Total 1 (1 unread)

1 D294C PATIENT2940 F 67Y Mary Back Health System 2005.04.12 09:42:47

SE: 1002
M: 1002
ACC# ACC2940
W 879 : L 575

Ready

start