

3308B

ROTATING TUBE STAND

The Universal 3308B rotating tube stand offers the versatility needed to facilitate a wide range of general purpose radiographic procedures. The sturdy tube mount assembly ensures stability and precision alignment.



Longitudinal lock, roll lock and vertical lock release button with positive indexing at 90 degrees for fast, accurate cross table positioning.

Key Features

- Pressure activated column rotation for fast, accurate positioning
- Electric vertical and longitudinal locks provide maximum holding power
- 180° column rotation and centering detent allow for rapid patient positioning
- 9in (23cm) of transverse travel for total table top accessibility



3308B

ROTATING TUBE STAND




Technical Specifications:

Ceiling height:	Minimum 96in (244cm) Maximum 116in (295cm)
Vertical travel:	Minimum 12in (30.5cm) Maximum 75in (190.5cm)
Longitudinal travel:	72in (183cm) with 8ft rail
Longitudinal rails length:	8ft standard
Transverse travel:	9in (23cm)
Column rotation:	$\pm 180^\circ$ with detents 0° , 90° , 90°
Tube rotation:	$+180^\circ$, -180°
ETL Listed	

Tube Stand Options:

- 140KHU, 200KHU or 300KHU X-ray tubes
- Manual, semi-automatic or fully automatic collimators
- Trunnion mount in place of standard platform mount
- 12ft, 10ft or 6ft longitudinal rails



DEL MEDICAL

Toll Free: (800) 800-6006
 Fax: (847) 288-7011
 Email: info@delmedical.com
 www.delmedical.com
 241 Covington Drive, Bloomingdale, IL 60108

Universal

Universal branded products are manufactured and distributed by UMG/Del Medical. Del Medical is an ISO 13485 & ISO 9001 Certified facility. Specifications are subject to change without notice.

UMG

Toll Free: (800) 261-9808
 Fax: (914) 835-6111
 Email: info@umgxray.com
 www.umgxray.com
 28 Calvert Street, Harrison, NY 10528

08.13 #9000-3308BUX